CeQur Simplicity™

Mealtime Insulin Patch

2U per click

IMPORTANT REMINDERS

- Please read all instructions before using your CeQur Simplicity™ Patch (or 'Patch').
- Use the Patch only with U-100 rapid-acting (mealtime) insulin: Humalog® or NovoLog®.
- Use your Patch for bolus insulin needs at mealtimes and snack times, or for correcting high blood glucose.
- Do not fill your Patch with long-acting insulin or any other insulin or medication.
- Do not directly expose your Patch to oils, lotions, or sunscreen.
- Do not wear your Patch for longer than 7 days.

Refer to User Guide for full instructions by scanning the QR code.

*Visit **myceqursimplicity.com** for more information.*

CeQur Simplicity™

Mealtime Insulin Patch

Dose Count Card

2U per click

Count by 2

"2...4...6...8"



Each click delivers 2 units of insulin

IMPORTANT

- The first click of buttons for a new patch delivers 0.3 units short of 2 units.
- If buttons lock during a dose, the last click of the buttons may deliver less than 2 units.
- If uncertain about your delivered dose, stop and check your blood glucose in 1 to 2 hours.

©2026 CeQur Corporation

AW-MTL-00872 Ver. 1