

# Count by 2

“2...4...6...8”



Each squeeze delivers  
2 units of insulin



## IMPORTANT

- The first squeeze of buttons for a new patch delivers 0.3 units short of 2 units.
- If buttons lock during a dose, the last squeeze of the buttons may deliver less than 2 units.
- If uncertain about your delivered dose, stop and check your blood glucose in 1 to 2 hours.



CeQur Simplicity™

4-Day Insulin Patch

# EXCITING NEWS!

Your patch is now FDA-cleared to  
be worn for up to **4 days**.



*Your monthly pack contains **8 patches**.*



CeQur Cares™  
[ We Care For You ]

Talk to your Healthcare Provider about your Prescription. If you have any questions, please call  
CeQur Cares™ at **1-888-552-3787**  
or visit us at **[www.myceqursimplicity.com](http://www.myceqursimplicity.com)**.

